

Purified anti-p62 (SQSTM1)

Mouse IgG1, kappa
O95A2

BioLegend
From Revvity



100 µg

Artikel-Nr.: B937202 | **BioLegend** | Herstellernummer: 937202

152,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICC

Clone: O95A2

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-p62 (SQSTM1)

Clone: O95A2

Isotype: Mouse IgG1, kappa

Other_Names: Sequestosome 1, p62, Phosphotyrosine-Independent Ligand For The Lck SH2 Domain Of 62 KDa, EBI3-Associated Protein Of 60 KDa, Ubiquitin-Binding Protein p62, Oxidative Stress Induced Like, Autophagy Receptor P62

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: ICC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=19428>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>