

## PE anti-p38 MAPK Phospho (Thr180/Tyr182)

Mouse IgG1, kappa  
A16016A

**BioLegend**  
From Revvity



100 tests

Artikel-Nr.: B690204 | **BioLegend** | Herstellernummer: 690204

**298,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** ICFC

**Clone:** A16016A

**Format:** PE

**Isotype:** kappa, Mouse IgG1

**Species of Origin:** Mouse

**Verpackung:** 100 Tests

PE anti-p38 MAPK Phospho (Thr180/Tyr182)

Clone: A16016A

Isotype: Mouse IgG1, kappa

Other\_Names: p38alpha, Mitogen-activated protein kinase 14, MAPK 14, CSBP

Size: 100 tests

Reactivity: Human, Mouse

Format: PE

App\_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=18747>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>