

PhosphoPair PLK-1 (Thr210) Antibody Set

3F8; 2A3

BioLegend
From Revvity



1 kit

Artikel-Nr.: B699956 | [BioLegend](#) | Herstellernummer: 699956

151,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: WB
Clone: 3F8; 2A3
Format: SET
Species of Origin: N/A
Verpackung: 1 Kit

PhosphoPair PLK-1 (Thr210) Antibody Set

Clone: 3F8; 2A3
Isotype:
Other_Names: Polo Like Kinase 1,STPK13, Cell Cycle Regulated Protein Kinase
Size: 1 kit

Reactivity: Human
Format: SET
App_Abbrev: WB
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1990>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>