

## TotalSeq-D1056 anti-human CD258 (LIGHT)

Mouse IgG2b, kappa  
T5-39

**BioLegend**  
From Revvity



10 µg

Artikel-Nr.: B318713 | BioLegend | Herstellernummer: 318713

**296,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** PG  
**Clone:** T5-39  
**Format:** TOTALSEQ\_D  
**Isotype:** kappa, Mouse IgG2b  
**Species of Origin:** Mouse  
**Verpackung:** 10 µg

TotalSeq-D1056 anti-human CD258 (LIGHT)

Clone: T5-39  
Isotype: Mouse IgG2b, kappa  
Other\_Names: LIGHT, TL5, TNFSF14, TR2 ligand, Ligand for herpesvirus entry mediator, tumor necrosis factor receptor-like-2 ligand, tumor necrosis factor ligand superfamily member 14 (TNFSF14)  
Size: 10 µg

Reactivity: Human  
Format: TOTALSEQ\_D  
App\_Abbrev: PG  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=21094>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>