

Purified anti-human CD222 (IGF2R) Recombinant

Mouse IgG1, kappa
QA19A18

BioLegend
From Revvity



100 µg

Artikel-Nr.: B364202 | BioLegend | Herstellernummer: 364202

310,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC, ICFC
Clone: QA19A18
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human CD222 (IGF2R) Recombinant

Clone: QA19A18
Isotype: Mouse IgG1, kappa
Other_Names: CD222, Insulin-like growth factor 2 Receptor (IGF2R), Cation-independent mannose 6-phosphate receptor (CIMPR), or Mannose-6 phosphate receptor (M6PR)
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: ICFC, FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20815>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>