

Purified anti-STAT4 Phospho (Tyr693)

Mouse IgG1, kappa
A19016A

BioLegend
From Revvity



100 µg

Artikel-Nr.: B941202 | BioLegend | Herstellernummer: 941202

264,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC, WB

Clone: A19016A

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-STAT4 Phospho (Tyr693)

Clone: A19016A

Isotype: Mouse IgG1, kappa

Other_Names: Signal Transducer and Activator Of Transcription 4, SLEB11, Signal Transducer and Activator of Transcription 4 Variant 3

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: WB, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20621>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>