

Purified anti-beta-Amyloid, 1-38 Recombinant

Rabbit IgG
QA19A28

BioLegend
From Revvity



100 µg

Artikel-Nr.: B943302 | BioLegend | Herstellernummer: 943302

259,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA, IHC-F, IHC-P, WB

Clone: QA19A28

Format: PURE

Isotype: Rabbit IgG

Species of Origin: Rabbit

Verpackung: 100 µg

Purified anti-beta-Amyloid, 1-38 Recombinant

Clone: QA19A28

Isotype: Rabbit IgG

Other_Names: APP, Amyloid Beta (A4) Precursor Protein, AD1, CVAP, Peptidase Nexin-II, Alzheimer Disease Amyloid Protein, Cerebral Vascular Amyloid Peptide, Protease Nexin-II, ABPP, APP-I, PN-II, AAA

Size: 100 µg

Reactivity:

Format: PURE

App_Abbrev: IHC-P, IHC-F, Direct ELISA, WB

Species of Origin: Rabbit

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20381>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>