

Ultra-LEAF Purified anti-human CXCL9 (MIG)

Mouse IgG1, kappa
A16023A

BioLegend
From Revvity



100 µg

Artikel-Nr.: B943503 | **BioLegend** | Herstellernummer: 943503

195,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Neut

Clone: A16023A

Format: ULEAF

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Ultra-LEAF Purified anti-human CXCL9 (MIG)

Clone: A16023A

Isotype: Mouse IgG1, kappa

Other_Names: Monokine induced by interferon gamma (MIG), Small inducible cytokine, subfamily B, member 9 (SCYB9)

Size: 100 µg

Reactivity: Human

Format: ULEAF

App_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20949>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>