

## Ultra-LEAF Purified anti-human CXCL9 (MIG)

Mouse IgG1, kappa  
A16023A

**BioLegend**  
From Revvity



1 mg

Artikel-Nr.: B943504 | BioLegend | Herstellernummer: 943504

**836,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** Neut  
**Clone:** A16023A  
**Format:** ULEAF  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 1 mg

Ultra-LEAF Purified anti-human CXCL9 (MIG)

Clone: A16023A  
Isotype: Mouse IgG1, kappa  
Other\_Names: Monokine induced by interferon gamma (MIG), Small inducible cytokine, subfamily B, member 9 (SCYB9)  
Size: 1 mg

Reactivity: Human  
Format: ULEAF  
App\_Abbrev: Neut  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20949>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>