

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Mouse IgG2a, kappa
A15129H

BioLegend
From Revvity



100 µg

Artikel-Nr.: B945403 | BioLegend | Herstellernummer: 945403

287,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Neut

Clone: A15129H

Format: ULEAF

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 100 µg

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Clone: A15129H

Isotype: Mouse IgG2a, kappa

Other_Names: Interferon-inducible T-cell alpha chemoattractant (ITAC), Small inducible cytokine, subfamily B, member 9B (SCYB9B),

Small inducible cytokine, subfamily B, member 11 (SCYB11)

Size: 100 µg

Reactivity: Human

Format: ULEAF

App_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20893>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>