

## Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Mouse IgG2a, kappa  
A15129H

**BioLegend**  
From Revvity



1 mg

Artikel-Nr.: B945404 | BioLegend | Herstellernummer: 945404

**964,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**App\_Abbrev:** Neut

**Clone:** A15129H

**Format:** ULEAF

**Isotype:** kappa, Mouse IgG2a

**Species of Origin:** Mouse

**Verpackung:** 1 mg

Ultra-LEAF Purified anti-human CXCL11 (ITAC)

Clone: A15129H

Isotype: Mouse IgG2a, kappa

Other\_Names: Interferon-inducible T-cell alpha chemoattractant (ITAC), Small inducible cytokine, subfamily B, member 9B (SCYB9B),

Small inducible cytokine, subfamily B, member 11 (SCYB11)

Size: 1 mg

Reactivity: Human

Format: ULEAF

App\_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=20893>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>