

Purified anti-human EMR3 Recombinant

Mouse IgG1, kappa
QA20A59

BioLegend
From Revvity



100 µg

Artikel-Nr.: B366102 | **BioLegend** | Herstellernummer: 366102

240,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: QA20A59
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-human EMR3 Recombinant

Clone: QA20A59
Isotype: Mouse IgG1, kappa
Other_Names: Adhesion G protein-coupled receptor E3; EGF module-containing mucin-like hormone receptor 3
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=21477>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>