

Ultra-LEAF Purified anti-mouse IL-4 Recombinant

Mouse IgG1, kappa
QA19A66

BioLegend
From Revvity



100 µg

Artikel-Nr.: B942803 | BioLegend | Herstellernummer: 942803

195,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Neut

Clone: QA19A66

Format: ULEAF

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Ultra-LEAF Purified anti-mouse IL-4 Recombinant

Clone: QA19A66

Isotype: Mouse IgG1, kappa

Other_Names: Interleukin-4, B-cell stimulating factor-1 (BSF-1), Ia inducing factor (IaIF), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 100 µg

Reactivity: Mouse

Format: ULEAF

App_Abbrev: Neut

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=21364>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>