

Purified anti-mouse/human Siglec-15

Mouse IgG1, kappa
MK4.5

BioLegend
From Revvity



100 µg

Artikel-Nr.: B164602 | **BioLegend** | Herstellernummer: 164602

165,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Direct ELISA, ELISA Capture, FC, ICC, ICFC

Clone: MK4.5

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-mouse/human Siglec-15

Clone: MK4.5

Isotype: Mouse IgG1, kappa

Other_Names: Sialic acid-binding Ig-like lectin 15, CD33 antigen-like 3

Size: 100 µg

Reactivity: Mouse, Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, ELISA Capture, Direct ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=21804>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>