

Purified anti-human CD198 (CCR8)

Rat IgG2b, kappa
S19016H

BioLegend
From Revvity



100 µg

Artikel-Nr.: B365802 | **BioLegend** | Herstellernummer: 365802

155,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: FC

Clone: S19016H

Format: PURE

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 100 µg

Purified anti-human CD198 (CCR8)

Clone: S19016H

Isotype: Rat IgG2b, kappa

Other_Names: CC-CKR-8, C-C CKR-8, CCR-8, CC chemokine receptor 8, CDw198, CDw198 antigen, CD198, chemokine (C-C motif) receptor 8, chemokine (C-C) receptor 8, TER1, ChemR1, and CKR-L1

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=21810>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>