



# LEGEND MAX Human Amyloid Beta (1-40) ELISA Kit



1 Pre-coated Plate

Artikel-Nr.: B449007 | [BioLegend](#) | Herstellernummer: 449007

**476,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** ELISA  
**Format:** LEGEND MAX  
**Species of Origin:** N/A  
**Verpackung:** 1 Pre-coated Plate

LEGEND MAX Human Amyloid Beta (1-40) ELISA Kit

Clone:

Isotype:

Other\_Names: AAA, ABETA, ABPP, AD1, APPI, CTFgamma, CVAP, PN-II, PN2, Amyloid beta A4 protein, preA4, protease nexin-II, peptidase nexin-II, beta-amyloid peptide, alzheimer disease amyloid protein, cerebral vascular amyloid peptide, APP, Amyloid precursor protein

Size: 1 Pre-coated Plate

Reactivity: Human  
Format: LEGEND MAX  
App\_Abbrev: ELISA  
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22285>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>