

## Purified anti-SKP1

Mouse IgG1, kappa  
P84E11C10

**BioLegend**  
From Revvity



25 µg

Artikel-Nr.: B605651 | [BioLegend](#) | Herstellernummer: 605651

**95,00 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

## Beschreibung

**App\_Abbrev:** ICC, WB  
**Clone:** P84E11C10  
**Format:** PURE  
**Isotype:** kappa, Mouse IgG1  
**Species of Origin:** Mouse  
**Verpackung:** 25 µg

Purified anti-SKP1

Clone: P84E11C10  
Isotype: Mouse IgG1, kappa  
Other\_Names: S-Phase Kinase Associated Protein 1, EMC19, OCP-II, OCP2, TCEB1L, p19a, P19skp1, SIII  
Size: 25 µg

Reactivity: Human, Mouse  
Format: PURE  
App\_Abbrev: WB, ICC  
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22169>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>