

PE anti-YAP1

Rat IgG2b, kappa
W19260A

BioLegend
From Revvity



100 tests

Artikel-Nr.: B603654 | **BioLegend** | Herstellernummer: 603654

243,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICFC

Clone: W19260A

Format: PE

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 100 Tests

PE anti-YAP1

Clone: W19260A

Isotype: Rat IgG2b, kappa

Other_Names: Yes1 Associated Transcriptional Regulator, YAP65, Yes-Associated Protein YAP65 Homolog, Transcriptional Coactivator YAP1, Yes Associated Protein 1, COB1, YAP2, YAP, YKI

Size: 100 tests

Reactivity: Human, Mouse

Format: PE

App_Abbrev: ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22448>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>