

TotalSeq-C0324 anti-mouse Hashtag 24

Rat IgG2a, kappa/Rat IgG2b, kappa
M1/42; 30-F11

BioLegend
From Revvity



10 µg

Artikel-Nr.: B113939 | BioLegend | Herstellernummer: 113939

310,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: PG
Clone: M1/42; 30-F11
Format: TOTALSEQ_C
Isotype: kappa, kappa/Rat IgG2b, Rat IgG2a
Species of Origin: Rat
Verpackung: 10 µg

TotalSeq-C0324 anti-mouse Hashtag 24

Clone: M1/42; 30-F11
Isotype: Rat IgG2a, kappa/Rat IgG2b, kappa
Other_Names: Mouse major histocompatibility complex H-2, MHC, T200, Ly-5, LCA
Size: 10 µg

Reactivity: Mouse
Format: TOTALSEQ_C
App_Abbrev: PG
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=23067>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>