

LEGEND MAX alpha-Synuclein Aggregate ELISA Kit

BioLegend
From Revvity



1 Pre-coated Plate

Artikel-Nr.: B449407 | **BioLegend** | Herstellernummer: 449407

488,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA
Format: LEGEND MAX
Species of Origin: N/A
Verpackung: 1 Pre-coated Plate

LEGEND MAX alpha-Synuclein Aggregate ELISA Kit

Clone:
Isotype:
Other_Names: alpha-synuclein oligomer, oligomeric alpha-synuclein, SNCA fibrillar aggregates
Size: 1 Pre-coated Plate

Reactivity: Human, Mouse
Format: LEGEND MAX
App_Abbrev: ELISA
Species of Origin: N/A

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22585>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>