

LEAF Purified anti-mouse SCF

Goat Polyclonal IgG

Poly5417

BioLegend
From Revvity



100 µg

Artikel-Nr.: B541703 | **BioLegend** | Herstellernummer: 541703

216,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: Block, WB

Clone: Poly5417

Format: LEAF

Isotype: Goat Polyclonal IgG

Species of Origin: Goat

Verpackung: 100 µg

LEAF Purified anti-mouse SCF

Clone: Poly5417

Isotype: Goat Polyclonal IgG

Other_Names: Kit Ligand, KITLG, Mast cell growth factor, MGF, Steel factor, SLF, Stem Cell Factor

Size: 100 µg

Reactivity: Mouse

Format: LEAF

App_Abbrev: Block, WB

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22807>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>