

Purified anti-alpha Smooth Muscle Actin

Mouse IgG2a, kappa
1A4

BioLegend
From Revvity



100 µg

Artikel-Nr.: B614852 | **BioLegend** | Herstellernummer: 614852

186,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICC, ICFC, IHC-P, WB

Clone: 1A4

Format: PURE

Isotype: kappa, Mouse IgG2a

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-alpha Smooth Muscle Actin

Clone: 1A4

Isotype: Mouse IgG2a, kappa

Other_Names: Actin, aortic smooth muscle (alphaSMA), alpha 2 actin (ACTA2), alpha-cardiac actin, growth-inhibiting gene 46, cell growth-inhibiting gene 46 protein

Size: 100 µg

Reactivity: Human, Mouse

Format: PURE

App_Abbrev: IHC-P, ICC, WB, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=23005>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>