


Caretta 96-Well PCR Platten, 2xR,G,G,B,V, farblos

PC Semi Skirt, High Profile, Eckschnitt A12

DNA-, DNase-, RNasefrei



Packung à 10 Stück

Artikel-Nr.: 711828 | Biozym

38,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Cyclertyp: Agilent-AriaMx qPCR System, Analytic Jena / Biometra-FlexCycler2, Analytic Jena / Biometra-qTOWERiris (touch), Analytic Jena / Biometra-T1, Analytic Jena / Biometra-TAdvanced 96, Analytic Jena / Biometra-TAdvanced 96G, Analytic Jena / Biometra-TAdvanced 96S, Analytic Jena / Biometra-TAdvanced 96SG, Analytic Jena / Biometra-TGradient, Analytic Jena / Biometra-TOne, Analytic Jena / Biometra-TProfessional Serie, Bio-Rad / MJ Research-C1000 (Touch), Bio-Rad / MJ Research-CFX96, Bio-Rad / MJ Research-CFX96 Touch, Bio-Rad / MJ Research-CFX96 Touch Deep Well, Bio-Rad / MJ Research-CFX Connect, Bio-Rad / MJ Research-DNA Engine Tetrad (2), Bio-Rad / MJ Research-DNA Engine, Bio-Rad / MJ Research-DNA Engine Dyad, Bio-Rad / MJ Research-DNA Engine Dyad Disciple, Bio-Rad / MJ Research-PC Tempo 96 / Deep Well, Bio-Rad / MJ Research-PTC 200, Bio-Rad / MJ Research-S1000, Bioer-GeneTouch 96, Bioer-GeneTouch Plus 96, Bioer-Lifeeco, Bioer-LifeTouch, Eppendorf-Mastercycler ep (gradient), Eppendorf-Mastercycler ep pro (S), Eppendorf-Mastercycler ep realplex, Eppendorf-Mastercycler nexus (gradient), Eppendorf-Mastercycler nexus GSX1, Eppendorf-Mastercycler nexus SX1, VWR-Collection UNO96

Aktionsangebote: Frühlingsaktion

Eckanschnitt: A12

Eigenschaften: rigid frame

Farbe: gemischt

Farbe Wells: farblos

Format: semi skirted

Kategorie: 96-Well Platte

Produktlinie: Caretta

Profil: High Profile

Verpackung: Packung à 10 Stück

Volumen: 0.2 ml