

APC anti-Nos2 (iNOS)

Rat IgG2b, kappa
W16030C

BioLegend
From Revvity



25 µg

Artikel-Nr.: B696807 | **BioLegend** | Herstellernummer: 696807

90,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: ICFC

Clone: W16030C

Format: APC

Isotype: kappa, Rat IgG2b

Species of Origin: Rat

Verpackung: 25 µg

APC anti-Nos2 (iNOS)

Clone: W16030C

Isotype: Rat IgG2b, kappa

Other_Names: Nitric oxide synthase 2, nitric oxide synthase 2A, NOS2A, HEP-NOS, inducible nitric oxide synthase (iNOS), Nitric oxide synthase, inducible, Inducible NOS

Size: 25 µg

Reactivity: Mouse

Format: APC

App_Abbrev: ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=22848>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>