

TotalSeq-B1286 anti-mouse CD284 (TLR4)

Rat IgG2a, kappa
SA15-21

BioLegend
From Revvity



10 µg

Artikel-Nr.: B145411 | **BioLegend** | Herstellernummer: 145411

222,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: PG
Clone: SA15-21
Format: TOTALSEQ_B
Isotype: kappa, Rat IgG2a
Species of Origin: Rat
Verpackung: 10 µg

TotalSeq-B1286 anti-mouse CD284 (TLR4)

Clone: SA15-21
Isotype: Rat IgG2a, kappa
Other_Names: Toll-like Receptor 4/MD2 complex, TLR4-MD2-complex, TLR4/MD2, CD284/MD2
Size: 10 µg

Reactivity: Mouse
Format: TOTALSEQ_B
App_Abbrev: PG
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=24243>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>