

Purified anti-human Granzyme A

Mouse IgG2b, kappa
W19350A

BioLegend
From Revvity



100 µg

Artikel-Nr.: B617152 | **BioLegend** | Herstellernummer: 617152

212,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA

Clone: W19350A

Format: PURE

Isotype: kappa, Mouse IgG2b

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human Granzyme A

Clone: W19350A

Isotype: Mouse IgG2b, kappa

Other_Names: Hanukah factor serine protease (HFSP), Cytolytic T cell-and natural killer cellspecific trypsin-like serine protease, Cytotoxic T-lymphocyte-associated serineesterase 3 (CTLA3)

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: Direct ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=23336>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>