

Purified anti-SLP76 Phospho (Tyr128)

Mouse IgG1, kappa
A20007B

BioLegend
From Revvity



100 µg

Artikel-Nr.: B618852 | **BioLegend** | Herstellernummer: 618852

308,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICC, ICFC

Clone: A20007B

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-SLP76 Phospho (Tyr128)

Clone: A20007B

Isotype: Mouse IgG1, kappa

Other_Names: Lymphocyte cytosolic protein 2 (LCP2), SH2 domain-containing leukocyte protein of 76 kDa (SLP76), IMD81

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: ICC, ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=23447>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>