

Purified anti-SOX17 Recombinant

Mouse IgG1, kappa
O92B6.Rec

BioLegend
From Revvity



100 µg

Artikel-Nr.: B628152 | **BioLegend** | Herstellernummer: 628152

285,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ICC, IHC-P, WB
Clone: O92B6.Rec
Format: PURE
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-SOX17 Recombinant

Clone: O92B6.Rec
Isotype: Mouse IgG1, kappa
Other_Names: Transcription factor SOX17, Transcription factor SOX 17, Transcription factor SOX-17, VUR3, SRY-Box17, SRY-Related HMG-Box Transcription Factor 17
Size: 100 µg

Reactivity: Human
Format: PURE
App_Abbrev: WB, ICC, IHC-P
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=24587>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>