

Caretta 96-Well PCR Platten, farblos, weiß

PC Full Skirt, Low Profile, Eckschnitt H1

DNA-, DNase-, RNasefrei



Packung à 10 Stück

Artikel-Nr.: 711809N | Biozym

51,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Cyclertyp: ABI-GeneAmp 2700, 2720, ABI-GeneAmp 9600, 9700, ABI-ProFlex, ABI-QuantStudio 3,5,6,7,12, ABI-SimpliAmp, ABI-Veriti 96-Well, Agilent-Mx3000P, Agilent-Mx3005P, Agilent-MX4000, Agilent-Surecycler 8800, Analytic Jena / Biometra-FlexCycler2, Analytic Jena / Biometra-qTOWER (G, touch), TOptical, Analytic Jena / Biometra-qTOWERiris (touch), Analytic Jena / Biometra-T1, Analytic Jena / Biometra-TAdvanced 96, Analytic Jena / Biometra-TAdvanced 96G, Analytic Jena / Biometra-TAdvanced 96S, Analytic Jena / Biometra-TAdvanced 96SG, Analytic Jena / Biometra-TGradient, Analytic Jena / Biometra-TOne, Analytic Jena / Biometra-TProfessional Serie, Bio-Rad / MJ Research-C1000 (Touch), Bio-Rad / MJ Research-CFX96 Touch Deep Well, Bio-Rad / MJ Research-DNA Engine Tetrad (2), Bio-Rad / MJ Research-DNA Engine, Bio-Rad / MJ Research-DNA Engine Dyad, Bio-Rad / MJ Research-DNA Engine Dyad Disciple, Bio-Rad / MJ Research-iCycler, iQ, iQ-5, Bio-Rad / MJ Research-MyiQ PCR System, Bio-Rad / MJ Research-PC Tempo 96 / Deep Well, Bio-Rad / MJ Research-PTC 200, Bio-Rad / MJ Research-S1000, Bio-Rad / MJ Research-T100, Bioer-GeneExplorer 96, Bioer-GeneTouch 96, Bioer-GeneTouch Plus 96, Bioer-Lifeeco, Bioer-LifeTouch, Eppendorf-Mastercycler ep (gradient), Eppendorf-Mastercycler ep pro (S), Eppendorf-Mastercycler ep realplex, Eppendorf-Mastercycler nexus (gradient), Eppendorf-Mastercycler nexus GSX1, Eppendorf-Mastercycler nexus SX1, peqLab-peqSTAR 96X, SensoQuest-LabCycler, VWR-Collection UNO96

Aktionsangebote: Frühlingsaktion

Eckanschnitt: H1

Eigenschaften: rigid frame

Farbe: Rahmen: farblos, weiß

Format: full skirted

Kategorie: 96-Well Platte

Produktlinie: Caretta

Profil: Low Profile

Verpackung: Packung à 10 Stück

Volumen: 0.1 ml