

FrameStar 96 PCR Plate, clear Wells, purple Frame

High Profile, cut H1

Non-Skirted



Packung à 50 Stück

Artikel-Nr.: AF4TI-0710 | Azenta | Herstellernummer: 4TI-0710

224,50 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Eckanschnitt: H1

Eigenschaften: rigid frame

Farbe: violett

Farbe Wells: farblos

Format: unskirted

Kategorie: 96-Well Platte

Produktlinie: FrameStar

Profil: High Profile

Verpackung: Packung à 50 Stück

Volumen: 0.2 ml

ABI® / LIFE TECHNOLOGIES / THERMO FISHER SCIENTIFIC

- GeneAmp® 2700 / 2720 / 9600 / 9700
- 7000, 7300, 7500, 7700, 7900 HAT
- 3100, 3130XL, 3500, 3500XL, 3730, 3730XL

AGILENT / STRATAGENE

- Surecycler 8800
- Robocycler Gradient
- Mx3000P,™ Mx3005P™
- Mx4000™

ANALYTIK JENA / BIOMETRA

- Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional (Standard/Basic) Gradient/XL
- Trobot 96, SpeedCycler2 (SP, SPR)
- qTOWER3 / G / touch, Toptica

BIOER TECHNOLOGIES

- Gene Touch 96

BIO-RAD

- C1000 Touch, S1000
- iCycler,™ MyCycler™, T100
- MyiQ,™ iCycler™ IQ / IQ 4 / IQ 5

BIO-RAD MJ RESEARCH

- Personal
- PTC100™ / 200™ / 220™ / 221™ / 225™ / 240™

Telefon: +49 (0) 51 52 / 90 20
Telefax: +49 (0) 51 52 / 20 70
E-Mail: support@biozym.com
Internet: www.biozym.com

Rechtliches
Registergericht: Amtsgericht Hannover
Registernummer: HRB 101682
UST-ID: DE 813739502

Biozym Scientific GmbH
Geschäftsführer
Dr. Sebastian Petri

- Chromo4™

CORBETT RESEARCH

-(Qiagen) Palm Cyclor

EPPENDORF

- MasterCycler® ep / ep gradient / Pro / Pro S / nexus / nexus gradient / nexus SX1 / nexus GSX1
- Mastercycler™ ep realplex

PEQLAB / VWR

- peqSTAR 96X

SENSOQUEST

- Labcycler

TAKARA

- Dice touch, Gradient

TECHNE

- Prime, PrimeG, Prime Elite, Prime Elite Satellite
- PCRmax Alpha cyclor 1 / 2 / 4
- TC412, TC512, Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene