


96 Well PCR Plate, clear

High Profile, cut H12

Non-Skirted



Packung à 50 Stück

Artikel-Nr.: AF4TI-0750 | Azenta | Herstellernummer: 4TI-0750

183,00 € *

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

Eckanschnitt: H12

Farbe: farblos

Format: unskirted

Kategorie: 96-Well Platte

Produktlinie: Azenta

Profil: High Profile

Verpackung: Packung à 50 Stück

Volumen: 0.2 ml

96-Well PCR Platte passend für:

ABI®/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- Thermal Cyclers (96 well standard block): Veriti | Proflex | Simplicamp | GeneAmp® 2700/2720/9600/9700
- qPCR Cyclers (96 well standard block): 7000 | 7300 | 7500 | 7700 | 7900 HT | QuantStudio™ 3/5/6/7/12K | ViiA7™

AGILENT/STRATAGENE

- Thermal Cyclers (96 well block): Surecycler 8800 | Robocycler Gradient
- qPCR Cyclers (96 well block): Mx3000P™ | Mx3005P™ | Mx4000™

ANALYTIK JENA/BIOMETRA

- Thermal Cyclers (96 well block): Flexcycler2 | T1 Thermocycler | Tgradient | Tone | Tadvanced | TProfessional (Standard/Basic) Gradient/XL | Trobot 96 | SpeedCycler2 (SP, SPR)
- qPCR Cyclers (96 well block): qTOWER3/G/touch | Toptical

BIOER TECHNOLOGIES

- Thermal Cyclers (96 well block): Gene Touch 96

BIO-RAD

Telefon: +49 (0) 51 52 / 90 20
Telefax: +49 (0) 51 52 / 20 70
E-Mail: support@biozym.com
Internet: www.biozym.com

Rechtliches
Registergericht: Amtsgericht Hannover
Registernummer: HRB 101682
UST-ID: DE 813739502

Biozym Scientific GmbH
Geschäftsführer
Dr. Sebastian Petri

- White plate ideal for qPCR, giving optimal signal-to-noise ratio for fluorescence-based assays

BIO-RAD MJ RESEARCH

- Thermal Cyclers (96 well block): PCT100™/200™/220™/221™/225™/240™
- qPCR Cyclers (96 well block): Chromo4™

CORBETT RESEARCH

- Thermal Cyclers (96 well block): (Qiagen) Palm Cycler

EPPENDORF

- Thermal Cyclers (96 well block): MasterCycler® ep/ep gradient/Pro/Pro S/nexus/nexus gradient/nexus SX1/nexus GSX1 | MasterCycler™ nexus X2/GX2/GX2e/X2e
- qPCR Cyclers (96 well block): MasterCycler™ ep realplex

PEQLAB/VWR

- Thermal Cyclers (96 well block): peqSTAR 96X

SENSOQUEST

- Thermal Cyclers (96 well block): Labcycler

TAKARA

- Thermal Cyclers (96 well block): Dice touch | Gradient

- **TECHNE**

- Thermal Cyclers (96 well block): Prime | PrimeG | Prime Elite | Prime Elite Satellite | PCRmax Alpha cycler 1/2/4 | TC412 | TC512 | Genius | Genius Quad | Touchgene | Touchgene Gradient | Flexigene