

## 384 Well PCR Plate, clear

cut A24  
Skirt



Packung à 50 Stück

Artikel-Nr.: AF4TI-1384 | Azenta | Herstellernummer: 4TI-1384

**196,50 € \***

\*zzgl. MwSt. [zzgl. Versandkosten](#)

### Beschreibung

**Eckanschnitt:** A24  
**Farbe:** farblos  
**Format:** full skirted  
**Kategorie:** 384-Well Platte  
**Produktlinie:** Azenta  
**Verpackung:** Packung à 50 Stück  
**Volumen:** 30 µl

## Diese 384-Well PCR Platte ist kompatibel mit folgenden Thermocyclern, qPCR Cyclern oder Sequencern

### ABI/LIFE TECHNOLOGIES/THERMO FISHER SCIENTIFIC

- Thermal Cyclers (384 well block): GeneAmp® 9700 | Verity | Proflex | Multiblock system
- qPCR Cyclers (384 well block): QuantStudio™ 5/6/7/12K | ViiA7™ | 7900 HT FAST
- Sequencers (384 well block): 3100 | 3130XL | 3500 | 3500XL | 3730 | 3730XL

### AGILENT/STRATAGENE

- Thermal Cyclers (384 well block): SureCycler 8800

### ANALYTIK JENA/BIOMETRA

- Thermal Cyclers (384 well block): Flexcycler2 | T1 Thermocycler | Tadvanced | TProfessional | Trobot 96
- qPCR Cyclers (384 well block): qTOWER3 84 / 84G

### BIOER TECHNOLOGIES

- Thermal Cyclers (384 well block): Gene Touch 384

### BIO-RAD

- Thermal Cyclers (384 well block): C1000 Touch | S1000

Telefon: +49 (0) 51 52 / 90 20  
Telefax: +49 (0) 51 52 / 20 70  
E-Mail: [support@biozym.com](mailto:support@biozym.com)  
Internet: [www.biozym.com](http://www.biozym.com)

Rechtliches  
Registergericht: Amtsgericht Hannover  
Registernummer: HRB 101682  
UST-ID: DE 813739502

Biozym Scientific GmbH  
Geschäftsführer  
Dr. Sebastian Petri

- qPCR Cyclers (384 well block): CFX384 Touch

#### **BIO-RAD MJ RESEARCH**

- Thermal Cyclers (384 well block): PCT200™/220™/221™/225™/240™

#### **CORBETT RESEARCH**

- Thermal Cyclers (384 well block): (Qiagen) Palm Cycler 384

#### **EPPENDORF**

- Thermal Cyclers (384 well block): MasterCycler® ep 384/Pro 384

#### **GE HEALTHCARE/AMERSHAM**

- Sequencers (384 well block): MegaBACE™ 4000

#### **PEQLAB/VWR**

- Thermal Cyclers (384 well block): peqSTAR 384X

#### **SENSOQUEST**

- Thermal Cyclers (384 well block): Labcycler

#### **TECHNE**

- Thermal Cyclers (384 well block): Prime | PrimeG | Prime Elite | Prime Elite Satellite | PCRmax Alpha cyclers 1/2/4 | TC412 | TC512 | Genius | Genius Quad | Flexigene