

Brilliant Violet 421 anti-human CD32B/C

Mouse IgG1, kappa
S18005H

BioLegend
From Revvity



100 tests

Artikel-Nr.: B398326 | **BioLegend** | Herstellernummer: 398326

264,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC
Clone: S18005H
Format: BV421
Isotype: kappa, Mouse IgG1
Species of Origin: Mouse
Verpackung: 100 Tests

Brilliant Violet 421 anti-human CD32B/C

Clone: S18005H
Isotype: Mouse IgG1, kappa
Other_Names: FCGR2B, FCgRIIB, FCGR2C, FCRIIC, IGFR2, CD32B, CD32C, FCG2, FCGR2, IGFR2, FcRII-c, CDW32
Size: 100 tests

Reactivity: Human
Format: BV421
App_Abbrev: FC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=25125>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>