

PE anti-NF-kappaB p65 Phospho (Ser529)

Mouse IgG1, kappa
A21012B

BioLegend
From Revvity



25 tests

Artikel-Nr.: B614153 | BioLegend | Herstellernummer: 614153

120,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ICFC

Clone: A21012B

Format: PE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 25 tests

PE anti-NF-kappaB p65 Phospho (Ser529)

Clone: A21012B

Isotype: Mouse IgG1, kappa

Other_Names: NF-kB p65, RelA,P65, NFKB3, Nuclear Factor NF-Kappa-B P65 subunit, Transcription Factor P65, NF-Kappa-B Transcription Factor P65, RELA Proto-Oncogene NF-KB subunit

Size: 25 tests

Reactivity: Human

Format: PE

App_Abbrev: ICFC

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=24616>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>