

Purified anti-human CCL7 (MCP-3)

Mouse IgG1, kappa
A15027A

BioLegend
From Revvity



100 µg

Artikel-Nr.: B630652 | [BioLegend](#) | Herstellernummer: 630652

332,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: Direct ELISA

Clone: A15027A

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 100 µg

Purified anti-human CCL7 (MCP-3)

Clone: A15027A

Isotype: Mouse IgG1, kappa

Other_Names: Monocyte chemotactic protein 3 (MCP-3), Small inducible cytokine subfamily A. member

Size: 100 µg

Reactivity: Human

Format: PURE

App_Abbrev: Direct ELISA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=24191>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>