

TotalSeq-B0387 anti-human TSLPR (TSLP-R)

Mouse IgG2a, lambda
1D3

BioLegend
From Revvity



10 µg

Artikel-Nr.: B322911 | **BioLegend** | Herstellernummer: 322911

296,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: PG
Clone: 1D3
Format: TOTALSEQ_B
Isotype: lambda, Mouse IgG2a
Species of Origin: Mouse
Verpackung: 10 µg

TotalSeq-B0387 anti-human TSLPR (TSLP-R)

Clone: 1D3
Isotype: Mouse IgG2a, lambda
Other_Names: Thymic stromal lymphopoietin protein receptor, cytokine receptor-like 2, CRL2, CRLF2, IL-XR, TSLP-R
Size: 10 µg

Reactivity: Human
Format: TOTALSEQ_B
App_Abbrev: PG
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=25386>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>