

2X RT-RPA Basic Master Mix, lyophilized



8 rxn of 50 µl

Artikel-Nr.: 351600300 | [biotechrabbit](#) | Herstellernummer: BC1600300

33,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Format: Lyophilisiert

Kategorie: Recombinase Polymerase Reaktion

Produkttyp: Kit

Verpackung: 8 rxn of 50 µl

Storage Temperature: 4°C or RT

Shipping Temperature: Ambient, Coolpack or Dry Ice

Minimal remaining shelf life (months): 6