

LYO-ready Exonuclease III, 100 U/μl



10000 U

Artikel-Nr.: 351500101 | [biotechrabbit](#) | Herstellernummer: BR1500101

72,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Format: Lyo-ready

Kategorie: Recombinase Polymerase Reaktion

Produkttyp: Komponente

Verpackung: 10000 U

Storage Temperature: -20°C

Shipping Temperature: Coolpack or Dry Ice

Minimal remaining shelf life (months): 6