

LYO-ready T4 UvsY Protein, 2 mg/ml



100 µg

Artikel-Nr.: 351500301 | [biotechrabbit](#) | Herstellernummer: BR1500301

170,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Format: Lyo-ready

Kategorie: Recombinase Polymerase Reaktion

Produkttyp: Komponente

Verpackung: 100 µg

Storage Temperature: -20°C

Shipping Temperature: Coolpack or Dry Ice

Minimal remaining shelf life (months): 6