

LYO-ready T4 Gene 32 Protein, 10 mg/ml



500 µg

Artikel-Nr.: 351500501 | [biotechrabbit](#) | Herstellernummer: BR1500501

288,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Format: Lyo-ready

Kategorie: Recombinase Polymerase Reaktion

Produkttyp: Komponente

Verpackung: 500 µg

Storage Temperature: -20°C

Shipping Temperature: Coolpack or Dry Ice

Minimal remaining shelf life (months): 6