

Purified anti-C/EBP beta

Mouse IgG2a, kappa
1H7

BioLegend
From Revvity



100 µg

Artikel-Nr.: B606202 | **BioLegend** | Herstellernummer: 606202

131,25 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion
App_Abbrev: ChIP, ICC, IP, WB
Clone: 1H7
Format: PURE
Isotype: kappa, Mouse IgG2a
Species of Origin: Mouse
Verpackung: 100 µg

Purified anti-C/EBP beta

Clone: 1H7
Isotype: Mouse IgG2a, kappa
Other_Names: CCAAT/enhancer-binding protein beta, Alpha-1-acid glycoprotein/enhancer binding protein, Interleukin 6-dependent DNA-binding protein (IL-6DBP), Nuclear factor of interleukin 6 (NF-IL6), Liver-enriched transcriptional activator protein (LAP)
Size: 100 µg

Reactivity: Mouse, Rat
Format: PURE
App_Abbrev: WB, ChIP, IP, ICC
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1984>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>