

Purified anti-human TNF-alpha

Mouse IgG1, kappa
MAb1

BioLegend
From Revvity



50 µg

Artikel-Nr.: B502801 | **BioLegend** | Herstellernummer: 502801

67,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: ELISA Capture, Neut, WB

Clone: MAb1

Format: PURE

Isotype: kappa, Mouse IgG1

Species of Origin: Mouse

Verpackung: 50 µg

Purified anti-human TNF-alpha

Clone: MAb1

Isotype: Mouse IgG1, kappa

Other_Names: Tumor necrosis factor-alpha, Cachectin, Necrosin, Macrophage cytotoxic factor (MCF), Differentiation inducing factor (DIF), TNFSF2

Size: 50 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture, Neut, WB

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1010>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>