

Purified anti-human IL-4

Rat IgG1, kappa
MP4-25D2

BioLegend
From Revvity



500 µg

Artikel-Nr.: B500802 | **BioLegend** | Herstellernummer: 500802

81,75 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

Aktionsangebote: Aktion

App_Abbrev: CyTOF®, ELISA, ICFC, IHC, Neut

Clone: MP4-25D2

Format: PURE

Isotype: kappa, Rat IgG1

Species of Origin: Rat

Verpackung: 500 µg

Purified anti-human IL-4

Clone: MP4-25D2

Isotype: Rat IgG1, kappa

Other_Names: Interleukin-4, Ia inducing factor (IaIF), B-cell stimulating factor-1 (BSF-1), Hodgkin's cell growth factor (HCGF), Mast cell growth factor-2 (MCGF-2), Macrophage fusion factor (MFF), T cell growth factor-2 (TCGF-2)

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA, CyTOF®, Neut, ICFC, IHC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=970>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>