

PE Goat anti-mouse IgG (minimal x-reactivity)

Goat Polyclonal IgG

Poly4053

BioLegend
From Revvity



200 µg

Artikel-Nr.: B405307 | **BioLegend** | Herstellernummer: 405307

176,00 € *

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FC

Clone: Poly4053

Format: PE

Isotype: Goat Polyclonal IgG

Species of Origin: Goat

Verpackung: 200 µg

PE Goat anti-mouse IgG (minimal x-reactivity)

Clone: Poly4053

Isotype: Goat Polyclonal IgG

Other_Names:

Size: 200 µg

Reactivity: Mouse

Format: PE

App_Abbrev: FC

Species of Origin: Goat

Weitere Informationen zum Produkt. Link zum Hersteller: <https://www.biolegend.com/Default.aspx?ID=6664&productid=1418>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>