

Purified anti-mouse Panendothelial Cell Antigen

Rat IgG2a, kappa
MECA-32



BioLegend
From Revvity

50 µg
Artikel-Nr.: B120501 | **BioLegend** | Hersteller-Nr.: 120501

57,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MECA-32
Format: PURE
Isotype: Rat IgG2a, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse Panendothelial Cell Antigen

Clone: MECA-32

Isotype: Rat IgG2a, kappa

Other_Names: Pan-endothelial Cell Antigen, MECA-32, Plasmalemma vesicle-associated protein

Size: 50 µg

Reactivity: Mouse

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=2967>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>