

Purified anti-human TNF-alpha

Mouse IgG1, kappa
MAb1



BioLegend
From Revvity

500 µg

Artikel-Nr.: B502802 | **BioLegend** | Hersteller-Nr.: 502802

188,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MAb1
Format: PURE
Isotype: Mouse IgG1, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human TNF-alpha

Clone: MAb1

Isotype: Mouse IgG1, kappa

Other_Names: Tumor necrosis factor-alpha, Cachectin, Necrosin, Macrophage cytotoxic factor (MCF), Differentiation inducing factor (DIF), TNFSF2

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=1010>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>