

Purified anti-human IL-8

Mouse IgG2a, kappa
H8A5



BioLegend
From Revvity

500 µg
Artikel-Nr.: B511502 | **BioLegend** | Hersteller-Nr.: 511502

226,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: ELISA Capture
Clone: H8A5
Format: PURE
Isotype: Mouse IgG2a, kappa
Species of Origin: Mouse
Verpackung: 500 µg

Purified anti-human IL-8

Clone: H8A5

Isotype: Mouse IgG2a, kappa

Other_Names: Interleukin-8, Neutrophil activating protein-1, Monocyte-derived neutrophil chemotactic factor (MDNCF), Neutrophil activating factor (NAF), Leukocyte adhesion inhibitor (LAI), Granulocyte chemotactic protein (GCP), CXCL8

Size: 500 µg

Reactivity: Human

Format: PURE

App_Abbrev: ELISA Capture

Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=4135>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>