

Alexa Fluor® 647 anti-mouse CD117 (c-kit)

Rat IgG2b, kappa
2B8



BioLegend
From Revvity

100 µg

Artikel-Nr.: B105818 | BioLegend | Hersteller-Nr.: 105818

184,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FC, ICFC
Clone: 2B8
Format: A647
Isotype: Rat IgG2b, kappa
Species of Origin: Rat
Verpackung: 100 µg

Alexa Fluor® 647 anti-mouse CD117 (c-kit)

Clone: 2B8

Isotype: Rat IgG2b, kappa

Other_Names: c-KIT, Stem Cell Factor Receptor (SCFR)

Size: 100 µg

Reactivity: Mouse

Format: A647

App_Abbrev: FC, ICFC

Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:

<https://www.biolegend.com/Default.aspx?ID=6664&productid=3132>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>