

Purified anti-mouse/human PNAd

Rat IgM, kappa
MECA-79



BioLegend
From Revvity

50 µg
Artikel-Nr.: B120801 | **BioLegend** | Hersteller-Nr.: 120801

57,00 €*

*zzgl. MwSt. zzgl. Versandkosten

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: MECA-79
Format: PURE
Isotype: Rat IgM, kappa
Species of Origin: Rat
Verpackung: 50 µg

Purified anti-mouse/human PNAd

Clone: MECA-79
Isotype: Rat IgM, kappa
Other_Names: PNAd, Peripheral Node Addressin, MECA-79
Size: 50 µg

Reactivity: Mouse, Human
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Rat

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2975>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>