

Purified anti-Cytochrome c

Mouse IgG2b, kappa
7H8.2C12



BioLegend
From Revvity

50 µL

Artikel-Nr.: B612503 | [BioLegend](#) | Hersteller-Nr.: 612503

85,00 €*

*zzgl. MwSt. [zzgl. Versandkosten](#)

Beschreibung

App_Abbrev: FA, FC, ICC, ICFC, IHC, IP, WB
Clone: 7H8.2C12
Format: PURE
Isotype: Mouse IgG2b, kappa
Species of Origin: Mouse
Verpackung: 50 µl

Purified anti-Cytochrome c

Clone: 7H8.2C12
Isotype: Mouse IgG2b, kappa
Other_Names: Cyt c
Size: 50 µL

Reactivity: Human, Mouse, Rat, Hamster
Format: PURE
App_Abbrev: FC, ICFC, ICC, IHC, IP, WB, FA
Species of Origin: Mouse

Weitere Informationen zum Produkt. Link zum Hersteller:
<https://www.biolegend.com/Default.aspx?ID=6664&productid=2002>

Unsere datenschutzrechtlichen Hinweise finden Sie unter: <https://www.biozym.com/datenschutzrechtliche-hinweise>